

LIFEPAK[®] 1000 DEFIBRILLATOR

Works like you work.™



- **Manual Override/ ECG Monitoring**
- **Easy to Use**
- **Longest ECG Viewing Time**
- **Powerful Battery**

TO CONTACT US:

Performance

The LIFEPAK 1000 defibrillator meets current guidelines:

- ADAPTIV™ Biphasic Technology: Escalating energy to 360J
- cprMAX™ Technology supports current guidelines
- CPR Feedback via CODE-STAT™ 7.0 data review software with Advanced CPR Analytics
- High Definition LCD Screen – 15cm diagonal
- 4 seconds of ECG on screen
- ECG monitoring with 3 lead cable
- High Capacity Battery of 4.5Amp/h – 440 200J shocks / 17 hrs of monitoring

Ruggedness

- Operating Temperature: 0°C to 50°C
- IP55: Excellent protection available against water and dust for an AED
- Robust: 1 meter drop test performed

Reliability

- Battery Gauge always visible, even when the AED is OFF
- Daily Self Tests
- Software allows you to configure the ID of your AED
- Motion detection during ECG analysis



For further information please contact your local Physio-Control representative or visit www.physio-control.com



HEADQUARTERS / MANUFACTURING

Physio-Control, Inc.
11811 Willows Road NE
P. O. Box 97006
Redmond, WA 98073-9706 USA
Tel. 425 867 4000
Fax. 425 867 4121
www.physio-control.com

SALES OFFICES

Physio-Control UK
Medtronic Ltd
Suite One, Sherbourne House
Croxley Business Park
Watford, Herts
WD18 8WW
Tel. 44 1923 212 213
Fax. 44 1923 241 004

Canada
Medtronic of Canada
6733 Kitimat Road
Mississauga, ON
L5N 1W3
Tel. 888 879 0977
Fax. 416 430 6115